

## Product datasheet for **TA357689**

### FAM161A Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	<b>Expected reactivity:</b> Cow, Dog, Horse, Human, Mouse, Pig, Rabbit, Rat <b>Homology:</b> Cow: 93%; Dog: 93%; Horse: 93%; Human: 100%; Mouse: 79%; Pig: 93%; Rabbit: 90%; Rat: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	65kDa
Gene Name:	family with sequence similarity 161 member A
Database Link:	<a href="#">NP_115556</a> <a href="#">Entrez Gene 84140 Human</a> <a href="#">Q3B820</a>
Background:	This gene belongs to the FAM161 family. It is expressed mainly in the retina. Mouse studies suggested that this gene is involved in development of retinal progenitors during embryogenesis, and that its activity is restricted to mature photoreceptors after birth. Mutations in this gene cause autosomal recessive retinitis pigmentosa-28.
Synonyms:	FLJ13305; MGC129982; MGC129983



[View online »](#)

## Product images:



WB Suggested Anti-FAM161A Antibody  
Titration: 1.0 ug/ml  
Positive Control: Placenta