

Product datasheet for **TA357693**

LSMD1 (NAA38) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-NAA38 antibody is: synthetic peptide directed towards the C-terminal region of Human LSMD1
Specificity:	Expected reactivity: Cow, Dog, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat, Zebrafish Homology: Cow: 100%; Dog: 100%; Guinea Pig: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Rabbit: 100%; Rat: 100%; Zebrafish: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	19 kDa
Gene Name:	N(alpha)-acetyltransferase 38, NatC auxiliary subunit
Database Link:	Entrez Gene 84316 Human Q9BRA0-2
Synonyms:	2010003I05Rik; AW214405; Lsm8



[View online »](#)

Product images:



Host: Rabbit
Target Name: LSMD1
Sample Tissue: Uterus Tumor Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: LSMD1
Sample Type: Uterus Tumor
Antibody Dilution: 1.0 μ g/ml