

## Product datasheet for **TA357702**


### **Tmem164 Rabbit Polyclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	<b>Expected reactivity:</b> Cow, Dog, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat <b>Homology:</b> Cow: 100%; Dog: 100%; Guinea Pig: 86%; Horse: 100%; Human: 100%; Mouse: 100%; Rabbit: 100%; Rat: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	33kDa
Gene Name:	transmembrane protein 164
Database Link:	<a href="#">NP_808260</a> <a href="#">Entrez Gene 209497 Mouse</a> <a href="#">Q6PHN7</a>
Background:	The function of this protein is unknown.
Synonyms:	bB360B22.3; FLJ20173; FLJ22679; OTTHUMP00000023848; RP13-360B22.2



[View online »](#)

**Product images:**

WB Suggested Anti-Tmem164 Antibody  
Titration: 1.0 ug/ml  
Positive Control: Mouse Kidney