

Product datasheet for **TA357711**

Bco2 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Dog, Guinea Pig, Horse, Human, Mouse, Rat Homology: Dog: 79%; Guinea Pig: 79%; Horse: 86%; Human: 100%; Mouse: 93%; Rat: 93%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	60kDa
Gene Name:	beta-carotene oxygenase 2
Database Link:	NP_573480 Entrez Gene 170752 Mouse Q99NF1
Background:	The function of this protein is unknown.
Synonyms:	B-DIOX-II; BCDO2; FLJ34464; OTTHUMP00000196941



[View online »](#)

Product images:



WB Suggested Anti-Bco2 Antibody
Titration: 1.0 ug/ml
Positive Control: Mouse Kidney