

Product datasheet for **TA357875**

Cks1b Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Cow, Dog, Goat, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat, Yeast, Zebrafish Homology: Cow: 100%; Dog: 100%; Goat: 79%; Guinea Pig: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Rabbit: 100%; Rat: 100%; Yeast: 79%; Zebrafish: 93%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	10kDa
Gene Name:	CDC28 protein kinase 1b
Database Link:	NP_058600 Entrez Gene 54124 Mouse P61025
Background:	Cks1b binds to the catalytic subunit of the cyclin dependent kinases and is essential for their biological function.
Synonyms:	CKS-1; CKS1; ckshs1; PNAS-16; PNAS-18; PNAS-143



[View online »](#)

Product images:



WB Suggested Anti-Cks1b Antibody
Titration: 1.0 ug/ml
Positive Control: Mouse Thymus