

## Product datasheet for **TA358037**

### LRRK1 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	<b>Expected reactivity:</b> Cow, Dog, Guinea Pig, Horse, Human, Mouse, Pig, Rabbit, Rat <b>Homology:</b> Cow: 92%; Dog: 79%; Guinea Pig: 86%; Horse: 100%; Human: 100%; Mouse: 100%; Pig: 93%; Rabbit: 100%; Rat: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	144kDa
Gene Name:	leucine-rich repeat kinase 1
Database Link:	<a href="#">NP_078928</a> <a href="#">Entrez Gene 79705 Human</a> <a href="#">Q38SD2</a>
Background:	The function of this protein remains unknown.
Synonyms:	FLJ23119; FLJ27465; KIAA1790; RIPK6; Roco1
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)

**Product images:**

WB Suggested Anti-LRRK1 Antibody  
Titration: 1.0 ug/ml  
Positive Control: HepG2 Whole Cell