

## Product datasheet for **TA358070**

### **FLJ43806 Rabbit Polyclonal Antibody**

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	<b>Expected reactivity:</b> Dog, Horse, Human, Mouse, Rabbit, Rat <b>Homology:</b> Dog: 93%; Horse: 100%; Human: 100%; Mouse: 93%; Rabbit: 86%; Rat: 93%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	24kDa
Database Link:	<a href="#">NP_963922</a> <a href="#">Entrez Gene 23254 Human</a> <a href="#">Q674X7</a>
Background:	The function of this protein remains unknown.

[View online »](#)

**Product images:**

WB Suggested Anti-FLJ43806 Antibody  
Titration: 1.0 ug/ml  
Positive Control: COLO205 Whole Cell