

Product datasheet for **TA358174**

C18orf1 (LDLRAD4) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Horse, Human, Mouse, Rabbit, Rat Homology: Horse: 86%; Human: 100%; Mouse: 93%; Rabbit: 86%; Rat: 93%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	26kDa
Gene Name:	low density lipoprotein receptor class A domain containing 4
Database Link:	NP_852148 Entrez Gene 753 Human O15165-4
Background:	The function of this protein remains unknown.
Synonyms:	8230401C20Rik; A430083H02; A430108L08Rik; C18orf1; D18Ert653e; D330030L18Rik
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

Product images:



WB Suggested Anti-C18orf1 Antibody
Titration: 1.0 ug/ml
Positive Control: U937 Whole Cell