

## Product datasheet for **TA358227**

### Mago nashi homolog 2 (MAGOHB) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	<b>Expected reactivity:</b> Cow, Dog, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat, Zebrafish <b>Homology:</b> Cow: 100%; Dog: 93%; Guinea Pig: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Rabbit: 100%; Rat: 100%; Zebrafish: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	16kDa
Gene Name:	mago homolog B, exon junction complex core component
Database Link:	<a href="#">NP_060518</a> <a href="#">Entrez Gene 55110 Human</a> <a href="#">Q96A72</a>
Background:	MAGOHB is involved in mRNA splicing and in the nonsense-mediated decay (NMD) pathway.
Synonyms:	FLJ10292; mago; magoh; MAGOH2; MGN2
Protein Pathways:	Spliceosome



[View online »](#)

**Product images:**



WB Suggested Anti-MAGOHB Antibody  
Titration: 1.0 ug/ml  
Positive Control: A549 Whole Cell