

Product datasheet for **TA358295**

TBRG4 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Cow, Human, Pig, Rat Homology: Cow: 86%; Human: 100%; Pig: 79%; Rat: 79%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	58kDa
Gene Name:	transforming growth factor beta regulator 4
Database Link:	NP_112162 Entrez Gene 9238 Human Q969Z0-2
Background:	TBRG4 may play a role in cell cycle progression.
Synonyms:	CPR2; FASTKD4; H_TD2522F11.8; KIAA0948; OTTHUMP00000159556



[View online »](#)

Product images:



WB Suggested Anti-TBRG4 Antibody
Titration: 1.0 ug/ml
Positive Control: MCF7 Whole Cell