

Product datasheet for **TA358346**


Sarcalumenin (SRL) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Cow, Dog, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat Homology: Cow: 83%; Dog: 92%; Guinea Pig: 93%; Horse: 83%; Human: 100%; Mouse: 86%; Rabbit: 93%; Rat: 92%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	54kDa
Gene Name:	sarcalumenin
Database Link:	NP_001092284 Entrez Gene 6345 Human Q86TD4-2
Background:	SRL may be involved in the regulation of calcium transport.
Synonyms:	sarcalumenin



[View online »](#)

Product images:

WB Suggested Anti-SRL Antibody
Titration: 1.0 ug/ml
Positive Control: MDA-MB-435S Whole Cell