

Product datasheet for **TA358443**

FYB1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Dog, Horse, Human, Mouse Homology: Dog: 100%; Horse: 79%; Human: 100%; Mouse: 75%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	96kDa
Gene Name:	FYN binding protein
Database Link:	NP_955367 Entrez Gene 2533 Human E9PBV9
Background:	The protein encoded by this gene is an adapter for the FYN protein and LCP2 signaling cascades in T-cells. The encoded protein is involved in platelet activation and controls the expression of interleukin-2. Three transcript variants encoding different isoforms have been found for this gene.
Synonyms:	ADAP; FYB-120/130; p120/p130; PRO0823; SLAP-130; SLAP130
Protein Families:	Druggable Genome



[View online »](#)

Product images:



WB Suggested Anti-FYB Antibody
Titration: 1.0 ug/ml
Positive Control: A549 Whole Cell