

Product datasheet for TA358483

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com

OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

EU: info-de@origene.com CN: techsupport@origene.cn

HIP1 Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: WB

Reactivity: Human Host: Rabbit

Clonality: Polyclonal

Immunogen: The immunogen is a synthetic peptide directed towards the C-terminal region of Human HIP1

Specificity: Expected reactivity: Cow, Dog, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat, Zebrafish

Homology: Cow: 100%; Dog: 100%; Guinea Pig: 100%; Horse: 100%; Human: 100%; Mouse:

100%; Rabbit: 100%; Rat: 100%; Zebrafish: 86%

Formulation: Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2%

sucrose.

Note that this product is shipped as lyophilized powder to China customers.

Concentration: lot specific

Purification: Affinity Purified
Conjugation: Unconjugated

Storage: For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small

aliquots to prevent freeze-thaw cycles.

Stability: Shelf life: one year from despatch.

Predicted Protein Size: 114kDa

Gene Name: huntingtin interacting protein 1

Database Link: NP 005329

Entrez Gene 3092 Human

O00291

Background: The product of this gene is a membrane-associated protein that functions in clathrin-

mediated endocytosis and protein trafficking within the cell. The encoded protein binds to the huntingtin protein in the brain; this interaction is lost in Huntington's disease. Alternative


splicing results in multiple transcript variants.

Synonyms: HIP-I; ILWEQ; MGC126506





Product images:



Host: Rabbit Target Name: HIP1 Sample Tissue: A549 Cell Lysate Antibody Dilution: 1.0µg/ml

Host: Rabbit Target Name: HIP1

Sample Type: A549 Whole Cell lysates

Antibody Dilution: 1.0ug/ml