

Product datasheet for **TA358516**

IL28B (IFNL3) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	This antibody is pan-IL28 specific, as there is 100% homology to the human IL28A sequence, per //www.uniprot.org/uniprot/Q8IZJ0 .

Expected reactivity: Dog, Horse, Human, Mouse, Pig, Rabbit, Rat

Homology: Dog: 93%; Horse: 86%; Human: 100%; Mouse: 79%; Pig: 86%; Rabbit: 85%; Rat: 77%

Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	22kDa
Gene Name:	interferon, lambda 3
Database Link:	NP_742151 Entrez Gene 282617 Human Q8IZJ9



[View online »](#)

Background:

This gene encodes a cytokine distantly related to type I interferons and the IL-10 family. This gene, interleukin 28A (IL28A), and interleukin 29 (IL29) are three closely related cytokine genes that form a cytokine gene cluster on a chromosomal region mapped to 19q13. Expression of the cytokines encoded by the three genes can be induced by viral infection. All three cytokines have been shown to interact with a heterodimeric class II cytokine receptor that consists of interleukin 10 receptor, beta (IL10RB) and interleukin 28 receptor, alpha (IL28RA).

Synonyms:

IL-28B; IL28B; IL28C

Protein Families:

Druggable Genome, Secreted Protein, Transmembrane

Protein Pathways:

Cytokine-cytokine receptor interaction, Jak-STAT signaling pathway

Product images: