

Product datasheet for TA358546


NPIP (NPIPA1) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Human Homology: Human: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2–8°C up to 1 week. For long term storage, store at –20°C in small aliquots to prevent freeze–thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	39kDa
Gene Name:	nuclear pore complex interacting protein family member A1
Database Link:	NP_008916 Entrez Gene 9284 Human Q9UND3
Background:	The function of this protein remains unknown.
Synonyms:	morpheus; NPIP



Product images:



WB Suggested Anti-NPIP Antibody
 Titration: 1.0 ug/ml
 Positive Control: Hela Whole Cell