

## Product datasheet for **TA358569**

### Gstt2 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of Mouse Gstt2
Specificity:	<b>Expected reactivity:</b> Dog, Human, Mouse, Rabbit, Rat <b>Homology:</b> Dog: 77%; Human: 100%; Mouse: 92%; Rabbit: 93%; Rat: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	28kDa
Gene Name:	glutathione S-transferase, theta 2
Database Link:	<a href="#">NP_034491</a> <a href="#">Entrez Gene 14872 Mouse</a> <a href="#">Q61133</a>
Background:	Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles.
Synonyms:	MGC182032



[View online »](#)

Product images:



Host: Rabbit  
Target Name: Gstt2  
Sample Type: Mouse Liver lysates  
Antibody Dilution: 1.0ug/ml