

Product datasheet for **TA358646**

COQ6 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Cow, Dog, Guinea Pig, Horse, Human, Mouse, Pig, Rabbit, Rat, Zebrafish Homology: Cow: 86%; Dog: 93%; Guinea Pig: 86%; Horse: 93%; Human: 100%; Mouse: 93%; Pig: 93%; Rabbit: 100%; Rat: 93%; Zebrafish: 77%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	48kDa
Gene Name:	coenzyme Q6, monooxygenase
Database Link:	NP_872286 Entrez Gene 51004 Human B7Z3K8
Background:	The function of this protein remains unknown.
Synonyms:	CGI-10
Protein Families:	Druggable Genome
Protein Pathways:	Metabolic pathways, Ubiquinone and other terpenoid-quinone biosynthesis



[View online »](#)

Product images:



WB Suggested Anti-COQ6 Antibody
Titration: 1.0 ug/ml
Positive Control: COLO205 Whole Cell