

Product datasheet for **TA358802**

STMN4 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the C terminal region of human STMN4
Specificity:	Expected reactivity: Human Homology: Cow: 100%; Dog: 100%; Guinea Pig: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Rat: 100%; Yeast: 90%; Zebrafish: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	25 kDa
Gene Name:	stathmin 4
Database Link:	NP_001269982.1 Entrez Gene 81551 Human Q9H169-2
Synonyms:	MGC111012; RB3



[View online »](#)

Product images:



Host: Rabbit
Target Name: STMN4
Sample Type: Nerve Fiber Tumor Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: STMN4
Sample Tissue: Human Shoulder Nerve Fiber
Tumor lysates
Antibody Dilution: 1 μ g/ml