

## Product datasheet for **TA358808**

### SPAG11A Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	<b>Expected reactivity:</b> Human <b>Homology:</b> Human: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	12kDa
Gene Name:	sperm associated antigen 11A
Database Link:	<a href="#">NP_001075021</a> <a href="#">Entrez Gene 653423 Human</a> <a href="#">Q6PDA7-2</a>
Background:	The function of this protein remains unknown.
Synonyms:	EDDM2A; EP2; HE2



[View online »](#)

Product images:



WB Suggested Anti-SPAG11A Antibody  
Titration: 1.0 ug/ml  
Positive Control: HCT15 Whole Cell