

Product datasheet for **TA358873**


TCL1B Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Human Homology: Human: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	15kDa
Gene Name:	T-cell leukemia/lymphoma 1B
Database Link:	NP_004909 Entrez Gene 9623 Human O95988
Background:	TCL1B enhances the phosphorylation and activation of AKT1 and AKT2.
Synonyms:	SYN-1; TCL1; TML1

[View online »](#)

Product images:



Host: Rabbit
 Target Name: TCL1B
 Antibody Dilution: 1.0ug/ml
 Sample Type: MCF7 cell lysate