

## Product datasheet for **TA358906**

### **MRFAP1 Rabbit Polyclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	<b>Expected reactivity:</b> Cow, Dog, Human, Pig, Rabbit <b>Homology:</b> Cow: 100%; Dog: 93%; Human: 100%; Pig: 100%; Rabbit: 85%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	15kDa
Gene Name:	Morf4 family associated protein 1
Database Link:	<a href="#">NP_150638</a> <a href="#">Entrez Gene 93621 Human</a> <a href="#">Q9Y605</a>
Background:	The function of this protein remains unknown.
Synonyms:	OTTHUMP00000155159; PAM14; PGR1



[View online »](#)

Product images:



Host: Rabbit  
Target Name: MRFAP1  
Antibody Dilution: 1.0ug/ml  
Sample Type: Jurkat cell lysate