

Product datasheet for **TA358938**


GGNBP1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the N-terminal region of Human GGNBP1
Specificity:	Expected reactivity: Horse, Human, Rabbit Homology: Horse: 86%; Human: 100%; Rabbit: 92%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	10kDa
Gene Name:	gametogenetin binding protein 1 (pseudogene)
Database Link:	AAR26430 Entrez Gene 449520 Human Q5YKI7
Background:	This gene is the ortholog of the mouse gametogenetin-binding protein 1 gene. In human, the open reading frame is disrupted by a nonsense mutation after 8-aa; consequently, this gene is currently considered to be a unitary pseudogene in human even though it is functional in other mammals.
Synonyms:	Ggnbp1



[View online »](#)

Product images:


WB Suggested Anti-GGNBP1 Antibody
 Titration: 1.0 ug/ml
 Positive Control: Fetal Liver