

Product datasheet for **TA358941**

PPM1H Rabbit Polyclonal Antibody


Product data:

| | |
|-------------------------|---|
| Product Type: | Primary Antibodies |
| Applications: | WB |
| Reactivity: | Human |
| Host: | Rabbit |
| Clonality: | Polyclonal |
| Immunogen: | The immunogen is a synthetic peptide directed towards the N-terminal region of Human PPM1H |
| Specificity: | Expected reactivity: Human |
| Formulation: | Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i> |
| Concentration: | lot specific |
| Purification: | Affinity Purified |
| Conjugation: | Unconjugated |
| Storage: | For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Stability: | Shelf life: one year from despatch. |
| Predicted Protein Size: | 56kDa |
| Gene Name: | protein phosphatase, Mg ²⁺ /Mn ²⁺ dependent 1H |
| Database Link: | NP_065751 Entrez Gene 57460 Human Q9ULR3 |
| Background: | The function of this protein remains unknown. |
| Synonyms: | ARHCL1; FLJ13253; KIAA1157; NERPP-2C; URCC2 |



[View online »](#)

Product images:



Host: Rabbit
Target Name: PPM1H
Sample Tissue: 293T Cell Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: PPM1H
Sample Type: 293T Whole Cell lysates
Antibody Dilution: 1.0ug/ml