

Product datasheet for **TA358948**

CDV3 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of Human CDV3
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	23kDa
Gene Name:	CDV3 homolog
Database Link:	Entrez Gene 55573 Human Q9UKY7-2
Background:	The function of this protein remains unknown.
Synonyms:	H41



[View online »](#)

Product images:



Host: Rabbit
Target Name: CDV3
Sample Tissue: Placenta Lysate
Antibody Dilution: 1.0µg/ml

Host: Rabbit
Target Name: CDV3
Sample Type: Placenta lysates
Antibody Dilution: 1.0ug/ml