

## Product datasheet for **TA359012**

### WDR43 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the N-terminal region of Human WDR43
Specificity:	<b>Expected reactivity:</b> Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	74kDa
Gene Name:	WD repeat domain 43
Database Link:	<a href="#">NP_055946</a> <a href="#">Entrez Gene 23160 Human</a> <a href="#">Q15061</a>
Background:	The function of this protein remains unknown.
Synonyms:	KIAA0007; NET12



[View online »](#)

## Product images:



Host: Rabbit  
 Target Name: WDR43  
 Sample Tissue: 721\_B Cell Lysate  
 Antibody Dilution: 1.0µg/ml

Host: Rabbit  
 Target Name: WDR43  
 Sample Type: 721\_B Whole Cell lysates  
 Antibody Dilution: 1.0ug/ml