

Product datasheet for **TA359028**


AGAP5 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of human AGAP5
Specificity:	Expected reactivity: Human Homology: Human: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	39kDa
Gene Name:	ArfGAP with GTPase domain, ankyrin repeat and PH domain 5
Database Link:	NP_001137472.1 Entrez Gene 729092 Human A6NIR3
Synonyms:	AGAP-5; CTGLF2; Em:AC073389.1



[View online »](#)

Product images:

Host: Rabbit
Target Name: CCNYL2
Sample Type: HepG2 Whole Cell lysates
Antibody Dilution: 1.0ug/ml