

Product datasheet for **TA359044**

C1orf183 (FAM212B) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of Human FAM212B
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	32kDa
Gene Name:	family with sequence similarity 212 member B
Database Link:	Entrez Gene 55924 Human Q9NTI7
Synonyms:	C1orf183



[View online »](#)

Product images:



Host: Rabbit
Target Name: FAM212B
Sample Tissue: Placenta Lysate
Antibody Dilution: 1.0µg/ml

Host: Rabbit
Target Name: FAM212B
Sample Type: Placenta lysates
Antibody Dilution: 1.0ug/ml