

## Product datasheet for **TA359074**


### KLHL33 Rabbit Polyclonal Antibody

#### Product data:

|                         |   |
|-------------------------|---|
| Product Type:           | Primary Antibodies  |
| Applications:           | WB  |
| Reactivity:             | Human   |
| Host:                   | Rabbit  |
| Clonality:              | Polyclonal  |
| Immunogen:              | The immunogen is a synthetic peptide directed towards the N-terminal region of Human KLHL33   |
| Specificity:            | <b>Expected reactivity:</b> Human   |
| Formulation:            | Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.                                       |
| Concentration:          | lot specific  |
| Purification:           | Affinity Purified   |
| Conjugation:            | Unconjugated  |
| Storage:                | For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Stability:              | Shelf life: one year from despatch.   |
| Predicted Protein Size: | 58kDa   |
| Gene Name:              | kelch like family member 33   |
| Database Link:          | <a href="#">NP_001103467</a><br><a href="#">Entrez Gene 123103 Human</a><br><a href="#">A6NCF5</a>                                      |
| Background:             | The function of this protein remains unknown.   |
| Synonyms:               | KLHL33  |



[View online »](#)

**Product images:**


Host: Rabbit  
 Target Name: KLHL33  
 Sample Tissue: MCF7 Cell Lysate  
 Antibody Dilution: 1.0µg/ml

Host: Rabbit  
 Target Name: KLHL33  
 Sample Type: MCF7 Whole Cell lysates  
 Antibody Dilution: 1.0ug/ml