

Product datasheet for **TA359121**

TBKBP1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the C-terminal region of Human TBKBP1
Specificity:	Expected reactivity: Cow, Dog, Horse, Human, Mouse, Pig, Rabbit, Rat Homology: Cow: 93%; Dog: 93%; Horse: 93%; Human: 100%; Mouse: 93%; Pig: 93%; Rabbit: 93%; Rat: 93%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	67kDa
Gene Name:	TBK1 binding protein 1
Database Link:	XP_005257919 Entrez Gene 9755 Human A7MCY6
Background:	TBKBP1 is an adaptor protein that binds to TBK1 and is part of the interaction network in the TNF/NFKB pathway.
Synonyms:	ProSAPiP2; SINTBAD



[View online »](#)

Product images:



Host: Rabbit
Target Name: TBKBP1
Sample Tissue: Jurkat Cell Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: TBKBP1
Sample Type: Jurkat Whole Cell lysates
Antibody Dilution: 1.0ug/ml