


Product datasheet for **TA359128**

TSGA10IP Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of Human TSGA10IP
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	44kDa
Gene Name:	testis specific 10 interacting protein
Database Link:	Entrez Gene 254187 Human Q3SY00-2
Background:	The function of this protein remains unknown.
Synonyms:	FAM161C; FLJ32880; MGC120465; MGC120466; MGC120468

[View online »](#)

Product images:


Host: Rabbit
 Target Name: TSGA10IP
 Sample Tissue: HT1080 Cell Lysate
 Antibody Dilution: 1.0µg/ml

Host: Rabbit
 Target Name: TSGA10IP
 Sample Type: HT1080 Whole Cell lysates
 Antibody Dilution: 1.0ug/ml