

## Product datasheet for **TA359141**

### C17orf107 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of Human C17orf107
Specificity:	<b>Expected reactivity:</b> Human <b>Homology:</b> Human: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	20kDa
Gene Name:	chromosome 17 open reading frame 107
Database Link:	<a href="#">NP_001139008</a> <a href="#">Entrez Gene 100130311 Human</a> <a href="#">Q6ZR85</a>
Background:	The function of this protein remains unknown.
Synonyms:	FLJ99427



[View online »](#)

## Product images:



Host: Rabbit  
Target Name: C17orf107  
Sample Tissue: Jurkat Cell Lysate  
Antibody Dilution: 1.0 $\mu$ g/ml

Host: Rabbit  
Target Name: C17ORF107  
Sample Tissue: Human Jurkat Whole Cell  
Antibody Dilution: 1 $\mu$ g/ml