

Product datasheet for **TA359286**

IAH1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the N-terminal region of human IAH1
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	17kDa
Gene Name:	isoamyl acetate-hydrolyzing esterase 1 homolog
Database Link:	NP_001034702.1 Entrez Gene 285148 Human Q2TAA2
Synonyms:	MGC102860



[View online »](#)

Product images:



Host: Rabbit
Target Name: IAH1
Sample Tissue: 721_B Cell Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: IAH1
Sample Type: 721_B Whole Cell lysates
Antibody Dilution: 1.0ug/ml