

Product datasheet for **TA359287**

ANKRD44 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of human ANKRD44
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	34kDa
Gene Name:	ankyrin repeat domain 44
Database Link:	NP_001182073.1 Entrez Gene 91526 Human Q8N8A2
Synonyms:	MGC21968; MGC70444; OTTHUMP00000205626; PP6-ARS-B



[View online »](#)

Product images:



Host: Rabbit
 Target Name: ANKRD44
 Sample Tissue: 721_B Cell Lysate
 Antibody Dilution: 1.0µg/ml

Host: Rabbit
 Target Name: ANKRD44
 Sample Type: 721_B Whole Cell lysates
 Antibody Dilution: 1.0ug/ml