

Product datasheet for **TA359430**

CRISP2 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the C-terminal region of Human CRIS2
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	26kDa
Gene Name:	cysteine rich secretory protein 2
Database Link:	XP_005249415 Entrez Gene 7180 Human P16562
Synonyms:	CRISP-2; CT36; GAPDL5; MGC111136; OTTHUMP00000016589; TPX1; TSP1



[View online »](#)

Product images:



Host: Rabbit
Target Name: CRISP2
Sample Tissue: NCI-H226 Cell Lysate
Antibody Dilution: 1.0µg/ml

Host: Rabbit
Target Name: CRISP2
Sample Type: NCI-H226 Whole Cell lysates
Antibody Dilution: 1.0ug/ml