

Product datasheet for **TA359431**

EIF5A Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the C-terminal region of Human IF5A1
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	20kDa
Gene Name:	eukaryotic translation initiation factor 5A
Database Link:	Entrez Gene 1984 Human P63241-2
Synonyms:	eIF-4D; EIF-5A; eIF-5A-1; eIF-5A1; EIF5A1; eIF5A1; MGC99547; MGC104255



[View online »](#)

Product images:



Host: Rabbit
Target Name: IF5A1
Sample Tissue: 721_B Cell Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: IF5A1
Sample Type: 721_B Whole Cell lysates
Antibody Dilution: 1.0 μ g/ml