

Product datasheet for TA359512

RTEL1 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the C-terminal region of Human RTEL1
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	109kDa
Gene Name:	regulator of telomere elongation helicase 1
Database Link:	Entrez Gene 51750 Human Q9NZ71-9
Background:	This gene encodes a DNA helicase which functions in the stability, protection and elongation of telomeres and interacts with proteins in the shelterin complex known to protect telomeres during DNA replication. Mutations in this gene have been associated with dyskeratosis congenita and Hoyerall-Hreidarsson syndrome. Read-through transcription of this gene into the neighboring downstream gene, which encodes tumor necrosis factor receptor superfamily, member 6b, generates a non-coding transcript. Alternative splicing results in multiple transcript variants encoding different isoforms.



[View online »](#)


Synonyms: bK3184A7.3; C20orf41; DKFZP434C013; KIAA1088; NHL; RTEL

Product images:



Host: Rabbit
 Target Name: RTEL1
 Sample Tissue: 721_B Cell Lysate
 Antibody Dilution: 1.0µg/ml

Host: Rabbit
 Target Name: RTEL1
 Sample Type: 721_B Whole Cell lysates
 Antibody Dilution: 1ug/ml



Host: Rabbit
 Target Name: RTEL1
 Sample Tissue: MCF7 Cell Lysate
 Antibody Dilution: 1.0µg/ml

Host: Rabbit
 Target Name: RTEL1
 Sample Type: MCF7 Whole Cell lysates
 Antibody Dilution: 1.0ug/ml