

Product datasheet for TA359542

SSH3 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the N-terminal region of Human SSH3
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	43 kDa
Gene Name:	slingshot protein phosphatase 3
Database Link:	Entrez Gene 54961 Human Q8TE77-4
Background:	The ADF (actin-depolymerizing factor)/cofilin family is composed of stimulus-responsive mediators of actin dynamics. ADF/cofilin proteins are inactivated by kinases such as LIM domain kinase-1. The SSH family appears to play a role in actin dynamics by reactivating ADF/cofilin proteins in vivo.
Synonyms:	FLJ10928; FLJ20515; FLJ42240; hSSH-3L; SSH-3L; SSH3L



[View online »](#)

Product images:



Host: Rabbit
Target Name: SSH3
Sample Tissue: COLO205 Cell Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: SSH3
Sample Type: COLO205 Whole Cell
Antibody Dilution: 1.0ug/ml