

Product datasheet for **TA359627**

DUSP13 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-DUSP13 antibody is: synthetic peptide directed towards the N-terminal region of Human DS13B
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	21 kDa
Gene Name:	dual specificity phosphatase 13
Database Link:	XP_005269947.1 Entrez Gene 51207 Human Q9UII6



[View online »](#)

Background:

Members of the protein-tyrosine phosphatase superfamily cooperate with protein kinases to regulate cell proliferation and differentiation. This superfamily is separated into two families based on the substrate that is dephosphorylated. One family, the dual specificity phosphatases (DSPs) acts on both phosphotyrosine and phosphoserine/threonine residues. This gene encodes different but related DSP proteins through the use of non-overlapping open reading frames, alternate splicing, and presumed different transcription promoters. Expression of the distinct proteins from this gene has been found to be tissue specific and the proteins may be involved in postnatal development of specific tissues. A protein encoded by the upstream ORF was found in skeletal muscle, whereas the encoded protein from the downstream ORF was found only in testis. In mouse, a similar pattern of expression was found. Multiple alternatively spliced transcript variants were described, but the full-length sequence of only some were determined.

Synonyms:

BEDP; DUSP13A; DUSP13B; FLJ32450; MDSP; OTTHUMP00000019875; SKRP4; TMDP

Product images:

Host: Rabbit
Target Name: DS13B
Sample Tissue: 721_B Cell Lysate
Antibody Dilution: 1.0µg/ml

WB Suggested Anti-DS13B antibody Titration: 1 ug/mL
Sample Type: Human 721_B Whole Cell