

## Product datasheet for **TA359776**

### NHLH2 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-NHLH2 antibody is: synthetic peptide directed towards the N-terminal of Human HEN2
Specificity:	<b>Expected reactivity:</b> Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	14 kDa
Gene Name:	nescient helix-loop-helix 2
Database Link:	<a href="#">XP_006710729.1</a> <a href="#">Entrez Gene 4808 Human</a> <a href="#">Q02577</a>
Synonyms:	bHLHa34; HEN2; KIAA0490; NSCL-2; NSCL2; OTTHUMP00000013748



[View online »](#)

Product images:



Host: Rabbit  
Target Name: HEN2  
Sample Tissue: ACHN Cell Lysate  
Antibody Dilution: 1.0µg/ml

Host: Rabbit  
Target Name: HEN2  
Sample Type: ACHN Whole Cell  
Antibody Dilution: 1.0ug/ml