

Product datasheet for TA359790

SALL2 Rabbit Polyclonal Antibody


Product data:

| | |
|-------------------------|---|
| Product Type: | Primary Antibodies |
| Applications: | WB |
| Reactivity: | Human |
| Host: | Rabbit |
| Clonality: | Polyclonal |
| Immunogen: | The immunogen is a synthetic peptide directed towards the N terminal region of human SALL2 |
| Specificity: | Expected reactivity: Human |
| Formulation: | Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i> |
| Concentration: | lot specific |
| Purification: | Affinity purified |
| Conjugation: | Unconjugated |
| Storage: | For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Stability: | Shelf life: one year from despatch. |
| Predicted Protein Size: | 110 kDa |
| Gene Name: | spalt like transcription factor 2 |
| Database Link: | NP_001278375.1 Entrez Gene 6297 Human Q9Y467 |
| Synonyms: | FLJ10414; FLJ55746; HSAL2; KIAA0360; p150(Sal2); SAL2; ZNF795 |
| Protein Families: | Druggable Genome |



[View online »](#)

Product images:



Host: Rabbit
Target Name: SALL2
Sample Type: A549 Cell Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: SALL2
Sample Tissue: A549 Whole Cell lysates
Antibody Dilution: 1 μ g/ml