

Product datasheet for **TA359801**

SNRP70 (SNRNP70) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the C-terminal region of Human RU17
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	48kDa
Gene Name:	small nuclear ribonucleoprotein U1 subunit 70
Database Link:	NP_003080 Entrez Gene 6625 Human P08621
Synonyms:	RNPU1Z; RPU1; SNRP70; U1-70K; U1AP; U1AP1; U1RNP; U170K



[View online »](#)

Product images:



Host: Rabbit
Target Name: RU17
Sample Tissue: Ovary Tumor Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: RU17
Sample Type: Ovary Tumor lysates
Antibody Dilution: 1.0ug/ml