

Product datasheet for **TA359819**


ZNF32 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the N-terminal region of Human ZNF32
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	30 kDa
Gene Name:	zinc finger protein 32
Database Link:	NP_008904.1 Entrez Gene 7580 Human P17041
Background:	ZNF32 may be involved in transcriptional regulation.
Synonyms:	C2H2-546; KOX30



[View online »](#)

Product images:

Host: Rabbit
Target Name: ZNF32
Sample Tissue: COLO205 Cell Lysate
Antibody Dilution: 1.0 μ g/ml

WB Suggested Anti-ZNF32 antibody Titration: 1
ug/mL
Sample Type: Human COLO205 Whole Cell