

Product datasheet for **TA359887**

SR140 (U2SURP) Rabbit Polyclonal Antibody


Product data:

| | |
|-------------------------|---|
| Product Type: | Primary Antibodies |
| Applications: | WB |
| Reactivity: | Human |
| Host: | Rabbit |
| Clonality: | Polyclonal |
| Immunogen: | The immunogen is a synthetic peptide directed towards the N-terminal region of Human SR140 |
| Specificity: | Expected reactivity: Human |
| Formulation: | Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i> |
| Concentration: | lot specific |
| Purification: | Affinity purified |
| Conjugation: | Unconjugated |
| Storage: | For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Stability: | Shelf life: one year from despatch. |
| Predicted Protein Size: | 113kDa |
| Gene Name: | U2 snRNP associated SURP domain containing |
| Database Link: | Entrez Gene 23350 Human O15042 |
| Synonyms: | 2610101N10Rik; AU023006; Sr140 |



[View online »](#)

Product images:



Host: Rabbit
Target Name: SR140
Sample Tissue: 721_B Cell Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: SR140
Sample Type: 721_B Whole Cell lysates
Antibody Dilution: 1.0ug/ml