

Product datasheet for **TA359972**

RPF1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the N-terminal region of Human RPF1
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	38kDa
Gene Name:	ribosome production factor 1 homolog
Database Link:	NP_079341 Entrez Gene 80135 Human Q9H9Y2
Background:	RPF1 may be required for ribosome biogenesis.
Synonyms:	BXDC5; RP11-118B23.1



[View online »](#)

Product images:



Host: Rabbit
Target Name: RPF1
Sample Tissue: Stomach Tumor Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: RPF1
Sample Type: Stomach Tumor lysates
Antibody Dilution: 1.0ug/ml