

## Product datasheet for **TA360198**

### PTBP1 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the N terminal region of human PTBP1
Specificity:	<b>Expected reactivity:</b> Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	57 kDa
Gene Name:	polypyrimidine tract binding protein 1
Database Link:	<a href="#">NP_002810.1</a> <a href="#">Entrez Gene 5725 Human</a> <a href="#">P26599</a>



[View online »](#)

**Background:**


This gene belongs to the subfamily of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs). The hnRNPs are RNA-binding proteins and they complex with heterogeneous nuclear RNA (hnRNA). These proteins are associated with pre-mRNAs in the nucleus and appear to influence pre-mRNA processing and other aspects of mRNA metabolism and transport. While all of the hnRNPs are present in the nucleus, some seem to shuttle between the nucleus and the cytoplasm. The hnRNP proteins have distinct nucleic acid binding properties. The protein encoded by this gene has four repeats of quasi-RNA recognition motif (RRM) domains that bind RNAs. This protein binds to the intronic polypyrimidine tracts that requires pre-mRNA splicing and acts via the protein degradation ubiquitin-proteasome pathway. It may also promote the binding of U2 snRNP to pre-mRNAs. This protein is localized in the nucleoplasm and it is also detected in the perinucleolar structure. Alternatively spliced transcript variants encoding different isoforms have been described.

**Synonyms:**

HNRNP-I; HNRNPI; HNRPI; MGC8461; MGC10830; pPTB; PTB; PTB-1; PTB-T; PTB2; PTB3; PTB4

**Protein Families:**

Druggable Genome

**Product images:**

Host: Rabbit  
Target Name: PTBP1  
Sample Tissue: HCT116 Cell Lysate  
Antibody Dilution: 1.0g/ml

Host: Rabbit  
Target Name: PTBP1  
Sample Tissue: Human HCT116 Whole Cell  
Antibody Dilution: 1.0ug/ml