

## Product datasheet for **TA360209**

### **RbAp48 (RBBP4) Rabbit Polyclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the C terminal region of human RBBP4
Specificity:	<b>Expected reactivity:</b> Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	48 kDa
Gene Name:	RB binding protein 4, chromatin remodeling factor
Database Link:	<a href="#">NP_001128727.1</a> <a href="#">Entrez Gene 5928 Human</a> <a href="#">Q09028</a>



[View online »](#)

**Background:**

This gene encodes a ubiquitously expressed nuclear protein which belongs to a highly conserved subfamily of WD-repeat proteins. It is present in protein complexes involved in histone acetylation and chromatin assembly. It is part of the Mi-2 complex which has been implicated in chromatin remodeling and transcriptional repression associated with histone deacetylation. This encoded protein is also part of co-repressor complexes, which is an integral component of transcriptional silencing. It is found among several cellular proteins that bind directly to retinoblastoma protein to regulate cell proliferation. This protein also seems to be involved in transcriptional repression of E2F-responsive genes. Three transcript variants encoding different isoforms have been found for this gene.

**Synonyms:**

NURF55; RBAP48; RBBP-4

**Protein Families:**

Druggable Genome, Transcription Factors

**Product images:**

Host: Rabbit  
Target Name: RBBP4  
Sample Tissue: 293T Cell Lysate  
Antibody Dilution: 1.0g/ml

Host: Rabbit  
Target Name: RBBP4  
Sample Tissue: Human 293T Whole Cell  
Antibody Dilution: 1.0ug/ml