

Product datasheet for **TA360334**

Lefty2 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the N terminal region of human LEFTY2
Specificity:	Expected reactivity: Mouse
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	40 kDa
Gene Name:	left-right determination factor 2
Database Link:	NP_796073.1 Entrez Gene 320202 Mouse P57785
Synonyms:	EBAF; LEFTA; LEFTYA; MGC46222; TGF-beta-4; TGFB4



[View online »](#)

Product images:



Host: Rabbit
Target Name: Lefty2
Sample Type: Mouse Kidney Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: Lefty2
Sample Tissue: Mouse Kidney lysates
Antibody Dilution: 1.0 μ g/ml