

Product datasheet for **TA360395**

RHBDD3 Rabbit Polyclonal Antibody


Product data:

| | |
|-------------------------|---|
| Product Type: | Primary Antibodies |
| Applications: | WB |
| Reactivity: | Human |
| Host: | Rabbit |
| Clonality: | Polyclonal |
| Immunogen: | The immunogen is a synthetic peptide directed towards the C terminal region of human RHBDD3 |
| Specificity: | Expected reactivity: Human |
| Formulation: | Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. |
| Concentration: | lot specific |
| Purification: | Affinity purified |
| Conjugation: | Unconjugated |
| Storage: | For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Stability: | Shelf life: one year from despatch. |
| Predicted Protein Size: | 40 kDa |
| Gene Name: | rhomboïd domain containing 3 |
| Database Link: | NP_036397.1 Entrez Gene 25807 Human Q9Y3P4 |
| Synonyms: | C22orf3; HS984G1A; PTAG |
| Protein Families: | Druggable Genome, Transmembrane |



[View online »](#)

Product images:



Host: Rabbit
Target Name: RHBDD3
Sample Tissue: Hela Cell Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: RHBDD3
Sample Tissue: Human Hela Whole Cell
Antibody Dilution: 1.0ug/ml