

Product datasheet for **TA360451**

SIDT1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the N terminal region of human SIDT1
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	94 kDa
Gene Name:	SID1 transmembrane family member 1
Database Link:	NP_001295279.1 Entrez Gene 54847 Human Q9NXL6
Synonyms:	B830021E24Rik; FLJ20174; SID-1; SID1



[View online »](#)

Product images:



Host: Rabbit
Target Name: SIDT1
Sample Type: NCI-H226 Cell Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: SIDT1
Sample Tissue: Human NCI-H226 Whole Cell
Antibody Dilution: 1.0 μ g/ml