

## Product datasheet for TA360452

### CCHCR1 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the N terminal region of human CCHCR1
Specificity:	<b>Expected reactivity:</b> Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	86 kDa
Gene Name:	coiled-coil alpha-helical rod protein 1
Database Link:	<a href="#">NP_001099033.1</a> <a href="#">Entrez Gene 54535 Human</a> <a href="#">Q8TD31</a>
Synonyms:	C6orf18; HCR; MGC126371; MGC126372; Pg8; SBP

[View online »](#)

## Product images:



Host: Rabbit  
Target Name: CCHCR1  
Sample Type: Uterus Tumor Lysate  
Antibody Dilution: 1.0µg/ml

Host: Rabbit  
Target Name: CCHCR1  
Sample Tissue: Uterus Tumor lysates  
Antibody Dilution: 1ug/ml